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Use the plates containing the His⁺ transformants and screen for the Mut⁺ and Mut^S phenotype as described below.

1. Using a sterile toothpick, pick one colony and streak or patch one His⁺ transformant in a regular pattern on both an MM plate and an MD plate, making sure to patch the MM plate first.
2. Use a new toothpick for each transformant and continue until 100 transformants have been patched (2-3 plates).
3. To differentiate Mut⁺ from Mut^S, make one patch for each of the controls (GS115/His⁺ Mut^S Albumin and GS115/His⁺ Mut⁺ β-gal) onto the MD and MM plates.
4. Incubate the plates at 30°C for 2 days.
5. After 2 days or longer at 30°C, score the plates. Mut⁺ transformants will grow well on both MD and MM plates. Mut^S transformants will grow well on MD plates, but show little or no growth on the MM plates.

We recommend purifying your His⁺ transformants to ensure isolation of a pure clonal isolates. You may do this before or after testing for the Mut phenotype.

Replica-Plating Procedure

This procedure gives a lower rate of misclassifications, but it increases the overall Mut⁺/MutS screening procedure by 2 days. You will need equipment to replica-plate.

1. Using sterile toothpicks, patch 100 His⁺ transformant on MD plates (2-3 plates). For controls, make one patch from each of the strains GS115/His⁺ MutS Albumin and GS115/His⁺ Mut⁺β⁻gal onto the MD plates.
2. Incubate the plates at 28-30°C for 2 days.
3. After 2 days, replica-plate the patches from the MD plates onto fresh MM and MD plates to screen for MutS transformants.
4. Incubate the replica plates at 28-30°C for 2 days.
5. After 2 days at 28-30°C, score the replica plates. Look for patches that grow normally on the MD replica plates but show little or no growth on the MM replica plates. Including His⁺ Mut⁺ and His⁺ MutS control patches on each plate will provide examples of Mut⁺ and MutS phenotypes.

Screening by Functional Assay

Some researchers have used a functional assay to directly screen for high expressing Pichia recombinant clones without first screening for MutS or Mut⁺ phenotypes. If you elect to screen directly for high-expressing recombinants, be sure to also check the Mut phenotype. This will help you optimize expression of your recombinant clone.

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2.7 BMGY 和 BMMY

Buffered Glycerol complex Medium, Buffered Methanol complex Medium (含有 1% yeast extract, 2% peptone, 100mM potassium phosphate pH 6.0, 1.34%YNB; 4X10⁻⁵% 生物素 1% glycerol or 0.5% methanol)

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You will need either BMGY/BMMY (buffered complex glycerol or methanol medium), BMG/BMM (buffered minimal glycerol or methanol medium) or MGY/MM (minimal glycerol or minimal methanol medium) for expression. **BMG, BMM, BMGY, and BMMY are usually used for the expression of secreted proteins,** particularly if pH is important for the activity of your protein.

Unlike MGY and MM, they are all buffered media. Because these media are buffered with phosphate buffer, a wide range of pH values may be used to optimize production of your protein. BMGY/BMMY contain yeast extract and peptone which may help stabilize secreted proteins and prevent or decrease proteolysis of secreted proteins. Inclusion of yeast extract and peptone act as a "mixed feed" allowing better growth and biomass accumulation.

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yeast extract	1%
peptone	2%
dextrose (glucose)	2%
sorbitol	1 M
+agar	2%

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Pichia yeast culture details

The growth temperature of *Pichia pastoris* is 28-30°C for liquid cultures, plates, and slants. Growth above 32°C during induction can be detrimental to protein expression and can even lead to cell death. Other important facts: Doubling time of log phase Mut+ or MutS *Pichia* in YPD is ~2 hours

- ï Mut+ and MutS strains do not differ in growth rates unless grown on methanol
- ï Doubling time of log phase Mut+ *Pichia* in methanol medium (MM) is 4-6 hours
- ï Doubling time of log phase MutS *Pichia* in MM is ~18 hours
- ï One OD600 = ~5 x 10⁷ cells/ml

Note: Growth characteristics may vary depending on the recombinant strain.

To store cells for weeks to months, use YPD medium or YPD agar slants (see page 55).

- Streak for single colonies of GS115, KM71, or a His+ transformant on YPD.
- Transfer one colony to a YPD stab and grow for 2 days at 30°C.
- The cells can be stored on YPD for several weeks at +4°C.

To store cells for months to years, store frozen at -80°C.

- Culture a single colony of GS115, KM71, or a His+ transformant overnight in YPD.
- Harvest the cells and suspend in YPD containing 15% glycerol at a final OD600 of 50-100 (approximately 2.5-5.0 x 10⁹ cells/ml).
- Cells are frozen in liquid nitrogen or a dry ice/ethanol bath and then stored at -80°C.

After long-term storage at +4°C or -80°C, we recommend checking the His+ transformants for correct genotype and viability. Streak on MM, MD or MGY plates before using again.

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0. make 100% (M/V)TCA:454 ml H₂O/Kg TCA, maintain in dark bottle at 4 °C, Be careful, use gloves!

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For PAGE-SDS, resuspend samples in a minimal volume of sample buffer. (The presence of some TCA can give a yellow colour as a consequence of the acidification of the sample buffer ; titrate with 1N NaOH or 1M TrisHCl pH8.5 to obtain the normal blue sample buffer colour.)

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- 1)、 The first is adding to the culture medium of amino-acid rich supplement such as peptone or casamino acid, which reduce product degradation possibly by acting as excess substrates for one or more problem proteases.
 - 2)、 The second strategy to minimize proteolytic degradation is to optimize the culture pH between pH 3.0 and pH7.0.
 - 3)、 The third strategy involves elevating the level of ammonium in the culture broth adequately, which is explained by the hypothesis that protease activity is induced under nitrogen starvation.
- The fourth strategy is let the culture temperature 28 deg instead of 30 deg.